

LALLA CARE : AN AGENTIC AI BABY MONITORING SYSTEM .

CHUNK TITLE: Overview of the Lalla Care Website

CHUNK CONTENT:

Lalla Care is a baby monitoring website that helps parents track their baby in real time through a live video feed. It offers features such as sleep tracking, activity history, environmental monitoring, alerts, and a built-in assistant to answer questions. The website is designed for non-technical users and provides a simple, organized layout with six main pages: Dashboard, Live Feed, Sleep Logs, Activity, Alerts, and Settings. A chatbot assistant is available on all pages to help answer user questions about either baby care or how to use the website.

CHUNK TITLE: Main Navigation Bar Overview

CHUNK CONTENT:

The main navigation bar appears at the very top of every page. It contains the Lalla Care logo on the left and five navigation tabs on the right. The tabs are Dashboard, Live Feed, Sleep Logs, Activity, and Settings. Clicking any tab switches to its corresponding page. The active page is marked with a highlighted underline. This navigation bar serves as the primary way to move between sections of the website.

CHUNK TITLE: Navigation Tabs - Layout and Contents

CHUNK CONTENT:

The navigation tabs appear horizontally across the top of all pages. The tabs are:

Dashboard: Opens the home page with camera feed and daily summaries.

Live Feed: Opens a dedicated view of the camera with enhanced controls.

Sleep Logs: Displays sleep patterns, statistics, and history.

Activity: Shows all recorded activities such as sleep, wake, feeding, crying, and adult visits.

Settings: Allows configuration of camera URL, alert cooldown, and detection preferences.

These tabs remain fixed in place as the user navigates.

CHUNK TITLE: Dashboard Page – General Overview

CHUNK CONTENT:

The Dashboard is the main landing page and central hub of the website. It provides a real-time

view of the baby's environment, daily summary cards, quick actions, audio monitoring, environmental conditions, and recent activity. The Dashboard automatically loads first when the user enters the site or when clicking the Dashboard tab at the top.

CHUNK TITLE: Dashboard – Last Update Indicator

CHUNK CONTENT:

The Last Update Indicator is located in the top right of the Dashboard page. It displays the time of the most recent data update along with a small pulsing indicator. This helps users confirm that the monitoring system is actively refreshing and connected. No interaction is required; the indicator updates automatically.

CHUNK TITLE: Dashboard – Live Camera Feed Area

CHUNK CONTENT:

The Live Camera Feed occupies a large area on the left side of the Dashboard. It includes a camera header showing the camera name and a label indicating the feed is live. When connected, the video appears automatically. If no connection exists, a placeholder message such as “Establishing Connection...” is shown. This part of the Dashboard is for viewing only and cannot be interacted with directly.

CHUNK TITLE: Dashboard – Daily Summary Cards

CHUNK CONTENT:

The Daily Summary Cards are located below the camera feed on the left side. They show key daily statistics including total sleep duration, number of feeding sessions, number of diaper changes, and total sleep change compared to previous days. Each card includes a value and a percentage increase or decrease, displayed using up or down arrows. The cards are informational and do not require user interaction.

CHUNK TITLE: Dashboard – Control Buttons Area

CHUNK CONTENT:

Just below the summary cards on the left side are three control buttons: Night Light, Lullaby, and Alerts. The Night Light button toggles the nursery night light. The Lullaby button plays

soothing music. The Alerts button navigates directly to the Alerts page. These buttons provide quick access to commonly used actions.

CHUNK TITLE: Dashboard – Audio Monitor Card

CHUNK CONTENT:

The Audio Monitor Card is located on the right side of the Dashboard. It displays animated sound bars that represent the audio levels detected in the nursery. This visual representation helps users observe sound activity even when not listening directly. The card is informative only and has no buttons.

CHUNK TITLE: Dashboard – Temperature and Humidity Card

CHUNK CONTENT:

This card appears on the right side below the Audio Monitor. It displays the current temperature in Fahrenheit and humidity percentage in the nursery. A circular "Listen In" button is present at the center, enabling the user to activate live audio listening. Small trend indicators show how temperature and humidity have changed over different time intervals such as 10, 12, and 20 hours.

CHUNK TITLE: Dashboard – Recent Activity Feed

CHUNK CONTENT:

The Recent Activity Feed is at the bottom right of the Dashboard. It shows a chronological list of recent events such as sleep sessions, wake periods, crying detection, feeding sessions, and adult presence. Each item includes a description, time, and duration. The feed scrolls if many events exist.

CHUNK TITLE: Live Feed Page Overview

CHUNK CONTENT:

The Live Feed page provides an expanded camera view with controls for zoom, brightness, night mode, and recording. It is accessed by clicking the Live Feed tab in the top navigation bar. This page is designed for users who want full control over viewing and adjusting the camera feed.

CHUNK TITLE: Live Feed – Page Header

CHUNK CONTENT:

At the top of the Live Feed page is the header section, displaying the title "Live Camera Feed," the camera name, and the last update indicator. This helps users confirm they are viewing the most recent feed.

CHUNK TITLE: Live Feed – Enhanced Video Feed Section

CHUNK CONTENT:

The main video feed is displayed on the left side. It includes a larger view than the Dashboard and provides several user controls below it. These include a zoom slider ranging from 1x to 3x, a brightness slider ranging from 0 to 100 percent, a night mode toggle, and a recording button that turns red when activated. This section is designed for viewing and adjusting the camera feed.

CHUNK TITLE: Live Feed – Quick Action Buttons

CHUNK CONTENT:

Below the main video controls are additional action buttons. These include a snapshot button for taking a picture, a two-way audio button for speaking through the monitor, a night light toggle, and a lullaby playback button. These buttons allow users to interact with the baby's environment without navigating back to the Dashboard.

CHUNK TITLE: Live Feed – Status Cards Sidebar

CHUNK CONTENT:

On the right side of the Live Feed page is a vertical sidebar containing status cards. These cards display the baby's state (such as sleeping or awake), whether an adult is present, whether crying is detected, and the connection status for camera, audio, and system updates. Color-coded indicators help users assess status at a glance.

CHUNK TITLE: Sleep Logs Page Overview

CHUNK CONTENT:

The Sleep Logs page is accessed through the top navigation. It provides detailed information on

sleep patterns and statistics. This page helps parents understand how long the baby has slept, how often the baby wakes up, and the overall sleep quality.

CHUNK TITLE: Sleep Logs – Time Period Selector

CHUNK CONTENT:

At the top right of the Sleep Logs page is a time period selector. It allows users to view sleep data for Today, Week, or Month. Selecting a time period updates the statistics and charts shown on the page.

CHUNK TITLE: Sleep Logs – Statistics Cards

CHUNK CONTENT:

Five statistics cards are displayed near the top of the Sleep Logs page. They show total sleep duration, average session length, number of wake-ups, longest session duration, and an overall sleep quality score. These cards provide a condensed summary of the baby's sleep for the selected time period.

CHUNK TITLE: Sleep Logs – Sleep Pattern Chart

CHUNK CONTENT:

A large chart on the left side of the Sleep Logs page displays sleep patterns over time. Sleep periods are highlighted, while awake periods appear differently. Hour markers help users interpret when sleep occurred throughout the day.

CHUNK TITLE: Sleep Logs – Sleep Quality Score Gauge

CHUNK CONTENT:

On the right side of the Sleep Logs page is a circular gauge showing the baby's overall sleep quality score. This visual meter provides an intuitive understanding of how restful the baby's sleep has been.

CHUNK TITLE: Sleep Logs – Recent Sleep Sessions List

CHUNK CONTENT:

At the bottom of the Sleep Logs page is a scrollable list of recent sleep sessions. Each entry

includes start and end times, duration, date, and a quality label such as excellent, good, or poor. This gives users a detailed record of each sleep period.

CHUNK TITLE: Activity Page Overview

CHUNK CONTENT:

The Activity page provides a complete list of all events recorded. It helps parents understand the baby's entire day, including sleep periods, wake events, crying, feeding sessions, and adult visits. The page also includes summaries and insights based on selected time periods.

CHUNK TITLE: Activity Page – Time Period Selector

CHUNK CONTENT:

A time period selector in the top right allows the user to view activity for Today, Week, or Month. Selecting a time period updates the timeline and summaries on the page.

CHUNK TITLE: Activity Page – Activity Filter Buttons

CHUNK CONTENT:

Below the header are several filter buttons that allow users to narrow the view to specific activity types such as sleep, wake, crying, feeding, or adult visits. Clicking a filter highlights it and updates the list below to show only the chosen activity type.

CHUNK TITLE: Activity Page – Activity Timeline

CHUNK CONTENT:

The main part of the Activity page is a vertical timeline in the center. A line runs down the left side, with circular icons marking each activity. To the right of each icon is a card with details such as description, timestamp, duration, and additional notes. Users can scroll to view older activities.

CHUNK TITLE: Activity Page – Activity Summary Cards

CHUNK CONTENT:

At the bottom of the Activity page are three summary cards. These show the most active hours

of the day, the distribution of different activity types, and quick insights about patterns such as sleep consistency or feeding routine.

CHUNK TITLE: Alerts Page Overview

CHUNK CONTENT:

The Alerts page provides a full history of system alerts. These alerts notify parents of important events such as crying detection, unusual movement, adult presence, or environmental issues. This page helps users review what has happened and clear alerts when needed.

CHUNK TITLE: Alerts Page – Page Header

CHUNK CONTENT:

The top of the Alerts page contains the title "Alert History," a brief description, and a "Clear All" button positioned on the top right. This allows the user to remove all past alerts at once.

CHUNK TITLE: Alerts Page – Alert List

CHUNK CONTENT:

The Alert List occupies the main portion of the page. Each alert is displayed as a card with a colored left border indicating severity: red for critical, orange for warnings, and blue for informational. Each alert includes a message, timestamp, and a severity label. Users can scroll to browse all alerts.

CHUNK TITLE: Alerts Page – Empty State Display

CHUNK CONTENT:

When no alerts exist, an empty state appears at the center of the page. It shows a large clipboard icon, a message stating that no alerts are found, and a note explaining that the alert history is currently empty. This helps the user understand that no issues are present.

CHUNK TITLE: Settings Page Overview

CHUNK CONTENT:

The Settings page is accessed via the top navigation bar and allows users to configure the

monitoring system. Users can enter their camera URL, adjust alert frequency, enable or disable various detection options, and start or stop monitoring.

CHUNK TITLE: Settings Page – Camera Configuration Card

CHUNK CONTENT:

This card appears on the left side of the Settings page. It contains a field where users can enter their camera's URL and a field to adjust the alert cooldown time. Alert cooldown determines how many seconds must pass between repeated alerts. This card is essential for connecting the camera feed and controlling alert frequency.

CHUNK TITLE: Settings Page – Detection Settings Card

CHUNK CONTENT:

Located on the right side of the Settings page, this card includes checkboxes for enabling or disabling Baby Detection, Adult Detection, and Crying Detection. Users can customize the types of monitoring they want active.

CHUNK TITLE: Settings Page – Monitoring Control Card

CHUNK CONTENT:

At the bottom of the Settings page is a wide card with three main buttons: Start Monitoring, Stop Monitoring, and Save Settings. Save Settings applies changes. Start Monitoring activates the monitoring system using the saved preferences. Stop Monitoring halts all monitoring processes. When active, the Start button changes to show that monitoring is running.

CHUNK TITLE: Chatbot Assistant Overview

CHUNK CONTENT:

The chatbot assistant appears on every page as a floating circular button located at the bottom-right corner. Clicking the button opens a chat window. The assistant provides guidance on using the website and answers general baby care questions. The chatbot is always accessible.

CHUNK TITLE: Chatbot – Floating Button Details

CHUNK CONTENT:

The floating button stays anchored to the bottom-right corner. When closed, it displays a circular icon. A notification indicator may appear when new information is available. Users click this button to open or close the chatbot window.

CHUNK TITLE: Chatbot – Chat Window Header

CHUNK CONTENT:

When opened, the chat window appears with a header at the top. The header contains an icon, the title "BabyBot Assistant," a subheading stating the bot is ready to help, and a close button on the right side. Clicking the close button hides the chat window.

CHUNK TITLE: Chatbot – Messages Area

CHUNK CONTENT:

The central portion of the chatbot window displays the message history. Messages from the assistant appear on the left, and messages from the user appear on the right. Each message includes a timestamp. New messages automatically scroll into view. When the assistant is typing, an indicator is displayed.

CHUNK TITLE: Chatbot – Quick Questions Buttons

CHUNK CONTENT:

At the beginning of a conversation, four quick question buttons appear at the bottom of the messages area. These buttons allow users to ask common questions, such as inquiries about sleep, crying, feeding, or temperature.

CHUNK TITLE: Chatbot – Input Area

CHUNK CONTENT:

The bottom section of the chatbot contains a text input box with a placeholder that prompts users to type their question. A send button positioned to the right transmits the message. Users can also press the Enter key to send questions.

CHUNK TITLE: Frequently Asked Questions – Camera Setup

CHUNK CONTENT:

Users often ask how to connect their camera. The answer is: open the Settings page, find the Camera Configuration card, enter the camera URL into the designated field, and click Save Settings.

CHUNK TITLE: Frequently Asked Questions – Start Monitoring

CHUNK CONTENT:

To begin monitoring, users must go to the Settings page, ensure that the correct camera URL and detection options are configured, click Save Settings, and then press the Start Monitoring button. The button changes appearance when monitoring is active.

CHUNK TITLE: Frequently Asked Questions – Alert Levels

CHUNK CONTENT:

Alerts have three levels. Informational alerts indicate general updates. Warning alerts indicate something that needs attention. Critical alerts indicate urgent situations that require immediate action. These can be reviewed on the Alerts page.

CHUNK TITLE: Frequently Asked Questions – Viewing Sleep Patterns

CHUNK CONTENT:

Sleep patterns can be viewed on the Sleep Logs page. Users can choose Today, Week, or Month using the time selector. Statistics cards and charts show patterns, quality, and session history.

CHUNK TITLE: Frequently Asked Questions – Controlling Nursery Features

CHUNK CONTENT:

Users can toggle the night light or play a lullaby from the Dashboard using the control buttons or from the Live Feed page under Quick Actions. These features help users soothe or check on the baby remotely.

CHUNK TITLE: Frequently Asked Questions – Daily Summary Meaning

CHUNK CONTENT:

The Daily Summary shows today's totals for sleep duration, feeding frequency, diaper changes, and overall sleep time. Percentage changes compare today's values to previous days.

CHUNK TITLE: Frequently Asked Questions – Viewing All Activities

CHUNK CONTENT:

To view the baby's full activity timeline, users should go to the Activity page. They can filter by activity type and select a time period such as Today, Week, or Month.

CHUNK TITLE: Frequently Asked Questions – Alert Cooldown Explanation

CHUNK CONTENT:

Alert cooldown is the number of seconds between repeated alerts of the same type. This prevents receiving repeated notifications too frequently. The cooldown setting is located in the Camera Configuration card under Settings.

CHUNK TITLE: Frequently Asked Questions – Disabling Detection Types

CHUNK CONTENT:

Users can disable Baby Detection, Adult Detection, or Crying Detection through the Detection Settings card on the Settings page by unchecking the corresponding boxes.

CHUNK TITLE: Frequently Asked Questions – Clearing All Alerts

CHUNK CONTENT:

Users can remove all alerts at once by opening the Alerts page and selecting the Clear All button at the top right.

CHUNK TITLE: Frequently Asked Questions – Audio Monitor Function

CHUNK CONTENT:

The Audio Monitor on the Dashboard displays moving sound bars that reflect live audio levels from the nursery. This helps users understand noise activity even when sound is muted.

CHUNK TITLE: Frequently Asked Questions – Using Listen In Feature

CHUNK CONTENT:

The Listen In feature is located inside the Temperature and Humidity card on the Dashboard. Pressing the circular Listen In button activates live audio streaming from the nursery.

CHUNK TITLE: Frequently Asked Questions – Capturing a Snapshot

CHUNK CONTENT:

On the Live Feed page, the Capture Snapshot button below the video feed allows users to take a still picture from the camera. The snapshot is captured instantly.

CHUNK TITLE: Frequently Asked Questions – Status Indicators Meaning

CHUNK CONTENT:

Color-coded dots and borders are used throughout the website. A green dot indicates something is active or connected. A red dot indicates a critical alert. A pulsing dot means the system is receiving live updates.

CHUNK TITLE: Color and Visual Indicator Guide

CHUNK CONTENT:

Green and teal tones indicate normal operation and good conditions. Blue is used for informational updates. Orange indicates warnings that may require attention soon. Red indicates critical alerts. Alert cards display a colored border to match the alert severity. Icons for activities and states appear on Timeline cards and status cards to help differentiate events.

CHUNK TITLE: Quick Navigation Instructions

CHUNK CONTENT:

To view the baby: open Dashboard or Live Feed.
To check sleep patterns: open Sleep Logs.
To review today's events: open Activity.
To view alerts: open Alerts.
To change settings or connect a camera: open Settings.
To ask questions: click the assistant button on the bottom-right corner of the screen.

CHUNK TITLE: Real-Time Update Behavior**CHUNK CONTENT:**

The website updates data automatically. The camera feed refreshes continuously. Status cards show live information. Alerts appear in real time as events occur. Activity logs also update automatically. The Last Update indicator confirms that fresh data is being received.

CHUNK TITLE: Getting Started Steps for New Users**CHUNK CONTENT:**

New users should begin by opening the Settings page, entering their camera URL, setting an alert cooldown value, selecting which detections to enable, clicking Save Settings, and then clicking Start Monitoring. After monitoring begins, the Dashboard shows live data.

CHUNK TITLE: Recommended Tips for Best Use**CHUNK CONTENT:**

For the best experience, users should check the Dashboard frequently for quick updates, use Sleep Logs to track long-term patterns, check Activity for daily details, enable only needed detections to avoid unnecessary alerts, set an appropriate cooldown to prevent alert spam, and rely on the chatbot for assistance.